


LED Panel (openPBR)

Step by step tutorial.


Introduction

In a photo bioreactor setup the photon source is probably the most expensive and complicated part, we decided to build in a RGBW light source with an intensity of at least 300 μ E to make light-dependent growth at high densities possible.




30 15 0 15 30 45 60 75 90 90

Technical data:



- 5 m LED strip with 72W
- RGB LED and white LED alternating
- Up to 300µE intensity



More details you can find under GitHub: Datasheets

Parts list:

- 5 m RGBW LED strip 79,90 Euro

https://www.led-konzept.de/RGBW-LED-Streifen-12V-72W-pro-5m-RGB-warmweiss-500cm-mit-Litzenanschluss

- 120W Power adapter 64,49 Euro

https://www.led-konzept.de/MeanWell-LED-Netzteil-HLG-120H-12B-12V-120W-IP67

- Aluminium Plate 57 x 175 x 2 mm ~ 2 Euro
- 4x N-Channel MOSFET 0,67 Euro

https://www.reichelt.de/mosfet-n-ch-60v-50a-110w-0-018r-to220-stp55nf06-p257486.html?r=1

- 4x Aluminium cool bloc 0,83 Euro

https://www.reichelt.de/kuehlkoerper-25-mm-alu-60-k-w-sot-32-to-220-sk-12-25-sa-32-p227995.html?&trstct=pos 0

Total:150,39

Step by step tutorial:

1. Step: Aluminum Plate


Take a generic 2 mm thick aluminum plate and cut it to a 175 x 527 mm. Drill 5 mm holes for the sensor LED's.








The position of these holes is crucial they have to sit exactly on the opposite to the Opt101 sensor of the sensor bar.



Step 2: LED stripes

- Stick the Led stripes with the tape of 50 cm in a distance of 5 mm.
- Solder the 5 pin contacts with cable to the next LED stripe.







Step 3: Electronics


This LED panel is controlled by 4 N-Channel MOSFET. One MOSFET per color (R,G,B,W). These transistors divide the 12V in to 255 steps in the Arduino code.

- Use 2 -3 mm screws to attach the cooling blocks to the MOSFET
- Connect G (Gate) with a 100 ohm resistor to your Arduino digital pin red pin 3, green pin 4, blue pin 5, white pin 6.
- Connect S (Source) to ground (-).
- Connect D (Drain) to each color cable of the LED.
- Connect the 12V (+) of the LED to the network adapter.










4. Calibration of the LED Panel

To determine the intensity of your LED panel us a µmol meter. Put the LED Panel in the casing of the Open PBR. Paint a point in a raster on the front wall, every 25 mm from top to bottom and every 80 mm. Measure the intensity at every point and plot them in to a graph. In the case of this LED panel the maximum intensity is 300 μ mol/s*m² e.g. (300 μ E).







The software to control the LED's you find in the openPBR_control software.